

Product components

Components	Component number	Size-1 25 mL	Size-2 100 mL
Nuclease-free H ₂ O	RM20214	25mL	100 mL

Product Description

Suitable for molecular biology applications requiring nuclease-free water, such as reagent preparation and enzymatic reactions. No toxic agents like DEPC were used in production, preventing inhibition of enzymatic reactions.

Applications

Protein analysis, cell analysis, RNA analysis, epigenetics, protein expression and purification, genome editing, glycobiology and proteomics, cloning, and synthetic biology.

Storage

-20°C